

Scissors



**RS-5600**

McPHERSON-VANNAS
straight, sharp
cutting edge 5mm
comb. tip width 0.3mm
overall length 3"

**RS-5601**

McPHERSON-VANNAS
curved, sharp
cutting edge 5mm
comb. tip width 0.3mm
overall length 3"

**RS-5602**

McPHERSON-VANNAS
straight, sharp
cutting edge 5mm
comb. tip width 0.2mm
overall length 3"

**RS-5603**

McPHERSON-VANNAS
curved, sharp
cutting edge 5mm
comb. tip width 0.2mm
overall length 3"

**RS-5604**

McPHERSON-VANNAS
straight
sharp needle points
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"

**RS-5605**

McPHERSON-VANNAS
curved
sharp needle points
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"

RS-5630 McPHERSON-VANNAS
straight
cutting edge 6mm
comb. tip width 0.3mm
overall length 4"



RS-5631 McPHERSON-VANNAS
curved
cutting edge 6mm
comb. tip width 0.3mm
overall length 4"



RS-5610 VANNAS
straight
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"




RS-5611 VANNAS
curved
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"



RS-5618 VANNAS
angled on edge
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"



RS-5619 VANNAS
angled on flat
cutting edge 3mm
comb. tip width 0.3mm
overall length 3"



**RS-5620**

VANNAS
straight
cutting edge 3mm
comb. tip width 0.2mm
overall length 3"

**RS-5621**

VANNAS
curved
cutting edge 3mm
comb. tip width 0.2mm
overall length 3"

**RS-5650**

CASTROVIEJO
straight, sharp
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"

**RS-5658**

CASTROVIEJO
angular, sharp
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"

**RS-5668**

CASTROVIEJO
angular, blunt
cutting edge 9mm
comb. tip width 0.3mm
overall length 3¹/₂"

**RS-5671**

CASTROVIEJO
curved, sharp
cutting edge 10mm
comb. tip width 0.33mm
overall length 3³/₄"

RS-5673 CASTROVIEJO
curved, blunt
cutting edge 10mm
comb. tip width 0.33mm
overall length 3¹/₂"



RS-5674 CASTROVIEJO
13mm blades
angled on flat, blunt
overall length 4¹/₂"



RS-5675 NOYES
slightly curved, sharp
cutting edge 13mm
comb. tip width 0.4mm
overall length 4¹/₂"




RS-5676 NOYES
straight, sharp
cutting edge 13mm
comb. tip width 0.4mm
overall length 4¹/₂"



RS-5677 NOYES
straight, sharp-blunt
overall length 4¹/₂"



RS-5678 NOYES
straight, blunt-blunt
overall length 4¹/₂"



RS-5680

Micro Dissecting
straight, 8mm blades
flat handle
comb. tip width 0.33mm
overall length 6"

**RS-5681**

Micro Dissecting
curved, 8mm blades
flat handle
comb. tip width 0.33mm
overall length 6"

**RS-5692**

Micro Dissecting
straight, 6mm blades
round handle
comb. tip width 0.33mm
overall length 6"



RS-5693

McCLURE
Micro Dissecting
curved, 6mm blades
round handle
comb. tip width 0.33mm
overall length 6"

**RS-5788**

NOYES
angled on flat
12mm blades
length 4^{3/4}"

**RS-5802**

DeWECKER
sharp-sharp
7mm blades
length 4^{1/4}"



1-800-424-2984
www.roboz.com



**RS-5820**

NOYES
sharp-sharp
7mm blade
length 6"



RS-5840 Shown

RS-5840**RS-5841**

BONN
15mm blades
length 3^{1/2}
straight, sharp
curved, sharp



RS-5850 Shown

RS-5850**RS-5851**

Micro Dissecting
14mm blades
flat shanks & large rings
length 3^{3/4}"
straight, sharp
curved, sharp

- Micro Dissecting
23mm blades
flat shanks & large rings
length 4"
- RS-5852** straight, sharp
RS-5853 curved, sharp



RS-5852 Shown


- RS-5859** WILMER
angled on flat
sharp
blade length 15mm
length 4"



- RS-5865** Micro Dissecting
curved
flat shanks
large finger rings
blade length 19mm
length 4"



1-800-424-2984
www.roboz.com



**RS-5869**

BIRO
 angled on flat
 sharp outside edges
 blade length 17mm
 length 4"

**RS-5871**

GRADLE
 slightly curved
 sharp points
 blade length 12mm
 length 4"


**RS-5880****RS-5881**

Micro Dissecting
 20mm blades
 length 3 1/2"
 straight, sharp
 curved, sharp

RS-5880 Shown



- Micro Dissecting
length 4"
RS-5882 straight, sharp
blade length 24mm
RS-5883 curved, sharp
blade length 24mm



RS-5882 Shown

- Micro Dissecting
blade length 23mm
length 3 1/2"
RS-5910 straight, sharp
RS-5911 curved, sharp



RS-5910 Shown

- Micro Dissecting
straight, sharp
blade length 24mm
length 4 1/2"
RS-5912 standard
RS-5912L left handed version
RS-5912SC SURE CUT
RS-5914 carbide blades
RS-5914SC SURE CUT, carbide



RS-5912 Shown

1-800-424-2984
www.roboz.com





RS-5913 Shown

Micro Dissecting
curved, sharp
blade length 24mm
length 4 1/2"
RS-5913 standard
RS-5913L left handed version
RS-5913SC SURE CUT
RS-5915 carbide blades
RS-5915SC SURE CUT, carbide



RS-5918 Micro Dissecting
angled, sharp
blade length 25mm
length 4 1/2"



RS-5921 Micro Dissecting
curved on side, sharp
blade length 22mm
length 4 1/2"



RS-5925 LA GRANGE
double curve
sharp serrated
blade length 18mm
length 4 1/4"

RS-5928

O'BRIEN
angular
blade length 14mm
length 3³/₄"

**RS-5930**


McGUIRE
right
blade length 15mm
length 4"

**RS-5940**

STEVENS
blade length 19mm
length 4¹/₄"


RS-5941

straight, sharp
curved, sharp



RS-5940 Shown

1-800-424-2984
www.roboz.com





RS-5942 Shown

RS-5942

STEVENS

blade length 19mm

length 4¹/₄"

straight, blunt

RS-5943

curved, blunt



RS-5960 Shown

RS-5960

KNAPP

blade length 27mm

length 4"

straight, sharp-blunt

RS-5963

curved, sharp-blunt



RS-5962 Shown

RS-5962

KNAPP

blade length 27mm

length 4"


straight, blunt-blunt

RS-5965

curved, blunt-blunt



RS-5980 Micro Dissecting
blade length 22mm
length 4"
straight, blunt
RS-5981 curved, blunt



RS-5980 Shown

RS-5982 Micro Dissecting
blade length 25mm
length 4 1/2"
straight, blunt, standard
RS-5982L left handed version
RS-5984 straight, blunt, carbide
blades



RS-5983 Micro Dissecting
blade length 25mm
length 4 1/2"
curved, blunt, standard
RS-5983L left handed version
RS-5985 curved, blunt, carbide
blades



RS-5990 Shown



**RS-5998**

Micro Dissecting
angled
probe points
blade length 24mm
length 4^{1/2}"



RS-6010 Shown

RS-6010

METZENBAUM
length 5"
straight, blunt
blade length 25mm

RS-6010SC

SURE CUT
straight, blunt
blade length 25mm

RS-6013


curved, blunt
blade length 27mm

RS-6013SC

SURE CUT
curved, blunt
blade length 27mm

**RS-6015**

RAGNELL
curved
slender blunt points
blade length 27mm
length 5"



RS-6031

Micro Dissecting
blade length 40mm
length 5"
light curve, blunt



RS-6031 Shown

RS-6043

POTTS-SMITH
blade length 13mm
length 7"
45° angle



RS-6043 Shown

1-800-424-2984
www.roboz.com



KELLY
blade length 29mm
length 6¹/₄"
straight, sharp
curved, sharp

RS-6050

RS-6051



RS-6050 Shown

RS-6058

KELLY
angled, sharp
blade length 29mm
length 6¹/₄"



RS-6700 Light Operating
blade length 30mm
length 4³/₄"
RS-6706 straight, sharp-blunt
curved, sharp-blunt



RS-6700 Shown

RS-6702 Light Operating
blade length 30mm
length 4³/₄"
RS-6703 straight, sharp-sharp
curved, sharp-sharp



RS-6702 Shown

RS-6705 Light Operating
curved, blunt
blade length 30mm
length 4³/₄"



RS-6708 Light Operating
angled, sharp
blade length 30mm
length 4³/₄"



BROPHY
blade length 40mm
length 5 1/2"
straight, sharp
curved, sharp

RS-6730**RS-6731**

RS-6730 Shown

Light Operating
blade length 40mm
length 5"
straight, sharp-blunt
curved, sharp-blunt

RS-6750**RS-6751**

RS-6750 Shown



- Light Operating
blade length 40mm
length 5"
RS-6752 straight, sharp-sharp
RS-6753 curved, sharp-sharp



RS-6752 Shown

- Light Operating
blade length 40mm
length 5"
RS-6754 straight, blunt-blunt
RS-6755 curved, blunt-blunt



RS-6754 Shown

- DEAVER
blade length 43mm
length 5 1/2"
RS-6760 straight, sharp-blunt
RS-6761 curved, sharp-blunt



RS-6760 Shown



DEAVER

blade length 43mm

length 5 1/2"

straight, sharp-sharp

curved, sharp-sharp

RS-6762**RS-6763**

RS-6762 Shown

RS-6780

PERMA-SET

straight, sharp-blunt

blade length 45mm

length 5 1/2"

**RS-6782**

PERMA-SET

straight, sharp-sharp

blade length 45mm

length 5 1/2"



Operating
sharp-blunt

RS-6800	straight, length 4 ¹ / ₂ "
RS-6806	straight, length 5"
RS-6812	straight, length 5 ¹ / ₂ "
RS-6812SC	SURE CUT, straight length 5 ¹ / ₂ "
RS-6818	straight, length 6"
RS-6824	straight, length 6 ¹ / ₂ "
RS-6830	straight, length 7 ¹ / ₂ "
RS-6832	straight, carbide length 5 ¹ / ₂ "
RS-6837	curved, length 4 ¹ / ₂ "
RS-6843	curved, length 5 ¹ / ₂ "
RS-6843SC	SURE CUT, curved length 5 ¹ / ₂ "
RS-6849	curved, length 6 ¹ / ₂ "



RS-6812 Shown

1-800-424-2984
www.roboz.com




Operating
sharp-sharp

RS-6802	straight, length 4 ¹ / ₂ "
RS-6808	straight, length 5"
RS-6814	straight, length 5 ¹ / ₂ "
RS-6820	straight, length 6"
RS-6826	straight, length 6 ¹ / ₂ "
RS-6834	straight, carbide length 5 ¹ / ₂ "
RS-6839	curved, length 4 ¹ / ₂ "
RS-6845	curved, length 5 ¹ / ₂ "
RS-6851	curved, length 6 ¹ / ₂ "



RS-6814 Shown

1-800-424-2984
www.roboz.com



Operating
blunt-blunt

RS-6804	straight, length 4 ¹ / ₂ "
RS-6810	straight, length 5"
RS-6816	straight, length 5 ¹ / ₂ "
RS-6822	straight, length 6"
RS-6828	straight, length 6 ¹ / ₂ "
RS-6841	curved, length 4 ¹ / ₂ "
RS-6847	curved, length 5 ¹ / ₂ "
RS-6853	curved, length 6 ¹ / ₂ "



RS-6816 Shown

1-800-424-2984
www.roboz.com



MAYO
straight

RS-6870 blade length 43mm
length 5¹/₂"

RS-6870SC SURE CUT
blade length 43mm
length 5¹/₂"

RS-6872 blade length 53mm
length 6³/₄"

RS-6872L left handed version

RS-6872SC SURE CUT
blade length 53mm
length 6³/₄"

RS-6874SC SURE CUT
blade length 56mm
length 9"


RS-6880 carbide blades
blade length 43mm
length 5¹/₂"

RS-6882 carbide blades
blade length 53mm
length 6³/₄"



RS-6870 Shown

1-800-424-2984
www.roboz.com



MAYO
curved

RS-6871 blade length 43mm
length 5¹/₂"

RS-6871SC SURE CUT
blade length 43mm
length 5¹/₂"

RS-6873 blade length 53mm
length 6³/₄"

RS-6875SC SURE CUT
blade length 63mm
length 9"

RS-6881 carbide blades
blade length 43mm
length 5¹/₂"

RS-6883 carbide blades
blade length 53mm
length 6³/₄"



RS-6873 Shown

1-800-424-2984
www.roboz.com



MAYO-STILLE
rounded blades

RS-6890 straight
blade length 42mm
length 5¹/₂"

RS-6892 straight
blade length 49mm
length 6³/₄"

RS-6891 curved
blade length 42mm
length 5¹/₂"

RS-6893 curved
blade length 49mm
length 6³/₄"



RS-6890 Shown

1-800-424-2984
www.roboz.com



MAYO-NOBLE
blade length 48mm
length 6¹/₂"

RS-6920

straight

RS-6921

curved



RS-6920 Shown

DOYEN
abdominal
blade length 67mm
length 7"

RS-6930

straight

RS-6931

curved



RS-6930 Shown



Veterinary
heavy pattern
sharp-blunt, serrated
blade length 57mm
length 6¹/₂"

RS-6940

straight

RS-6941

curved



RS-6940 Shown

Veterinary
heavy pattern
blunt-blunt, serrated
blade length 57mm
length 6¹/₂"

RS-6942

straight

RS-6943


curved



RS-6942 Shown




RS-6947 DuBOIS
decapitation
curved, serrated
blade length 57mm
length 10^{1/2}"



50% of Actual Size

METZENBAUM
straight

RS-6950 blade length 39mm
length 5^{1/2}" (Lahey)
RS-6952 blade length 35mm
length 7" (Lahey)
RS-6952L left handed version



RS-6950 Shown

1-800-424-2984
www.roboz.com



METZENBAUM
extra delicate
curved
blade length 32mm

RS-6956


carbide blades
length 7"

RS-6957

length 7"

RS-6959

length 6"



RS-6959 Shown

1-800-424-2984
www.roboz.com




METZENBAUM
curved

RS-6951	blade length 39mm length 5 1/2" (Lahey)
RS-6951SC	SURE CUT blade length 39mm length 5 1/2" (Lahey)
RS-6953	blade length 35mm length 7" (Lahey)
RS-6953L	blade length 35mm length 7" (Lahey) Left Hand Version
RS-6953SC	SURE CUT blade length 35mm length 7" (Lahey)
RS-6954SC	SURE CUT blade length length 7"
RS-6955	carbide blades blade length 35mm length 7" (Lahey)
RS-6955SC	SURE CUT carbide blades blade length 35mm length 7" (Lahey)
RS-6973	blade length 49mm length 9" (Nelson)
RS-6958	carbide blades blade length 39mm length 5 1/2" (Lahey)
RS-6974	carbide blades blade length 49mm length 9" (Nelson)



RS-6951 Shown



RS-7010**RS-7011**

KAHN
blade length 32mm
straight
length 5³/₄"
curved
length 5³/₄"



RS-7010 Shown

RS-7030**RS-7031**

SIMS
blade length 48mm
length 8"
straight, sharp-sharp
curved, sharp-sharp

RS-7030 Shown
75% of Actual Size

RS-7041

DEAN
blades angled on flat
one serrated blade
slightly curved handles
blade length 21mm
length 6^{3/4}"

**RS-7045**

REYNOLDS
curved
tenotomy tips
blade length 35mm
length 6"



RS-7045 Shown



RS-7047


MARTIN
cartilage
serrated
blade length 48mm
length 7^{3/4"}



75% of Actual Size

RS-7055

HERCULES
heavy duty
serrated
length 7^{1/2"}



90% of Actual Size



RS-7060

Enterotomy
hook blade
length 8¹/₂"



75% of Actual Size

RS-7084

LISTER
bandage
length 4¹/₂"

RS-7085

length 5¹/₂"

RS-7084 Shown

1-800-424-2984
www.roboz.com



RS-7097

Extra Heavy Duty
hardened steel
serrated
length 7³/₄"



90% of Actual Size



RS-7098 ESMARCH
bandage
and plaster shears
length 8"



75% of Actual Size

RS-7099L All Purpose
bandage/utility
serrated, autoclavable
black plastic handles
length 7¹/₂"



90% of Actual Size



RS-7099S

All Purpose
bandage/utility
serrated, autoclavable
black plastic handles
length 5 1/2"



1-800-424-2984
www.roboz.com

