

Clips and Clamps



RS-5420

Micro Clip
straight
60gr pressure
clip width 0.8mm
jaw length 5mm

**RS-5422**

Micro Clip
straight
85gr pressure
clip width 0.75mm
jaw length 6mm

**RS-5424**

Micro Clip
straight
100gr pressure
clip width 1mm
jaw length 8mm

**RS-5426**

Micro Clip
straight
125gr pressure
clip width 1.5mm
jaw length 10mm

**RS-5431**

Micro Clip
slight curve
85gr pressure
clip width 1mm
jaw length 6mm

**RS-5433**

Micro Clip
strong curve
85gr pressure
clip width 1mm
jaw length 6mm

**RS-5435**

Micro Clip
45° angle
85gr pressure
clip width 1mm
jaw length 6mm



RS-5410

Clip Applying Forceps
for micro clips
length 5¹/₂"



RS-5440

Micro Clip
applying forceps
with lock
length 5³/₄"




1-800-424-2984
www.roboz.com



RS-5450

SCHWARTZ
clip applying forceps
for temporary clips
length 8"



75% of actual size

RS-5452

SCHWARTZ
temporary clip
straight, smooth
jaw width 1.7mm
jaw length 8mm
clip length 1"




RS-5457

SCHWARTZ
temporary clip
slight angle, smooth
jaw width 1.7mm
jaw length 8mm
clip length 1"



RS-5459

SCHWARTZ
temporary clip
strong angle, smooth
jaw width 1.7mm
jaw length 8mm
clip length 1"



RS-6470 Micro Vascular Clip
vessel diameter 1mm
occlusion pressure 10-15gr
length 8mm



RS-6472 Micro Vascular Clip
vessel diameter 1-2mm
occlusion pressure 15-20gr
length 11mm



RS-6474 Micro Vascular Clip
vessel diameter 2-3mm
occlusion pressure 20-30gr
length 14mm



RS-6488 Micro Vascular Approximator
with moveable clips
vessel diameter 1mm
occlusion pressure 10-15gr
length 8mm



RS-6490 Micro Vascular Approximator
with moveable clips
vessel diameter 1-2mm
occlusion pressure 15-20gr
length 11mm



RS-6492 Micro Vascular Approximator
with moveable clips
vessel diameter 2-3mm
occlusion pressure 20-30gr
length 14mm



RS-6496 Micro Clip Setting Forceps
length 4^{3/4"}



**RS-5470**

DIEFFENBACH
serrefine
straight, very delicate
length 1 ¹/₈"

**RS-5471**


DIEFFENBACH
serrefine
curved, very delicate
length 1 ¹/₈"

**RS-7438**

DIETRICH
micro bulldog clamps
straight
jaw length 8mm
jaw width 1.2mm
length 2"

**RS-7439**

DIETRICH
micro bulldog clamp
curved
jaw length 8mm
jaw width 1.2mm
length 2"



RS-7440	JOHNS HOPKINS bulldog clamp straight length 1 1/2"	
RS-7441	JOHNS HOPKINS bulldog clamp curved length 1 1/2"	
RS-7442	JOHNS HOPKINS bulldog clamp straight length 2"	
RS-7443	JOHNS HOPKINS bulldog clamp curved length 2"	
RS-7444	JOHNS HOPKINS bulldog clamp straight length 2 1/2"	
RS-7445	JOHNS HOPKINS bulldog clamp curved length 2 1/2"	
RS-7446	JOHNS HOPKINS bulldog clamp straight length 3"	
RS-7447	JOHNS HOPKINS bulldog clamp curved length 3"	




**RS-7500**

DeBAKEY
atraumatic bulldog clamp
straight
jaw length 2cm
overall length 7.5cm

**RS-7511**

DeBAKEY
atraumatic bulldog clamp
curved
jaw length 2cm
overall length 7.5cm



RS-7530

GLOVER
atraumatic bulldog clamp
straight
tension adjusting screw
jaw length 2cm
overall length 7cm



RS-7541

GLOVER
atraumatic bulldog clamp
curved
with tension adjusting screw
jaw length 2cm
overall length 7cm




1-800-424-2984
www.roboz.com



**RS-7760**

JONES
towel clamp
cross action
light weight
length 2¹/₄"



RS-7784 Shown

RS-7780**RS-7784**

BACKHAUS
towel clamp
length 3¹/₂"
length 5¹/₄"

1-800-424-2984
www.roboz.com

