

Wound Closures



**RS-7981**

Eye Needles
3/8 circle, cutting edge
size: 1, 2, 3, or 4
dozen

Please specify size when ordering

**RS-7982**

Eye Needles
1/2 circle, cutting edge
size: 1, 2, 3, or 4
sold by dozen

Please specify size when ordering

**RS-7983**


Eye Needles
1/2 circle, taper point
size: 1, 2, 3, or 4
sold by dozen

Please specify size when ordering

**RS-7985**

Regular Surgeon's
3/8 circle, cutting edge
sold by dozen

Please specify size when ordering




RS-7986 Regular Surgeon's
1/2 circle, cutting edge
sold by dozen



Please specify size when ordering

RS-7987 Fine Intestinal Needle
straight, taperpoint
sizes: 1 1/2"
sold by dozen




NON-STERILE

Cat.#	Size	Material		Length	Needle	Packaged
SUT-15-1	5-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-2	4-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-3	3-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool
SUT-15-4	2-0	Silk, Black, Braided	Non-Absorbable	100yd	None	Spool



STERILE


Cat.#	Size	Material		Length	Needle	Packaged
SUT-1234-13	6-0	Silk, Black, Braided	Non-Absorbable	18"	HE-6	Box/12
SUT-1228-13	6-0	Silk, Black, Braided	Non-Absorbable	18"	HE-7	Box/12
SUT-1073-11	6-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1073-21	5-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1074-21	5-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1091-23	5-0	Silk, Black, Braided	Non-Absorbable	18"	P-12	Box/12
SUT-570-31	4-0	Gut, Chromic	Absorbable	18"	C-12	Box/36
SUT-566-31	4-0	Gut, Chromic	Absorbable	18"	C-13	Box/36
SUT-1066-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-22	Box/36
SUT-1073-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-1	Box/36
SUT-1074-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1076-31	4-0	Silk, Black, Braided	Non-Absorbable	18"	C-14	Box/36
SUT-1203-33	4-0	Silk, Black, Braided	Non-Absorbable	18"	P-13	Box/12
SUT-1123-31	4-0	Silk, Black, Braided	Non-Absorbable	30"	V-20	Box/36
SUT-566-41	3-0	Gut, Chromic	Absorbable	18"	C-13	Box/36
SUT-1074-41	3-0	Silk, Black, Braided	Non-Absorbable	18"	C-13	Box/36
SUT-1076-41	3-0	Silk, Black, Braided	Non-Absorbable	18"	C-14	Box/36
SUT-1123-41	3-0	Silk, Black, Braided	Non-Absorbable	30"	V-20	Box/36

1-800-424-2984
www.roboz.com



NEEDLE CHART FOR SUTURES WITH NEEDLES

REVERSE CUTTING
3/8 CIRCLEREVERSE CUTTING
1/2 CIRCLETAPER
1/2 CIRCLE

**RS-7069**

Wire Cutting
serrated
angled with notch
length 5"

**RS-5950**

SPENCER
stitch, very delicate
length 3 1/2"

**RS-7074****RS-7076**

LITTAUER
stitch
length 4 1/2"
length 5 1/2"

RS-7074 Shown

**RS-6400**

KALT
with spring catch
length 5 1/2"



RS-6403

PATTON
straight without lock, serrated
length 4 $\frac{1}{2}$ "
serrated



RS-6410

CASTROVIEJO
length 5 $\frac{3}{4}$ "
straight without catch
RS-6412
straight with catch



RS-6412 Shown

RS-6416

RS-6417

CASTROVIEJO
extra delicate
0.6mm wide jaws
to hold 7-0 to 10-0
sutures
length 5 $\frac{3}{4}$ "
straight with catch
RS-6417
curved with catch



RS-6417 Shown





RS-6421 Shown

RS-6421

McPHERSON
length 4"
curved
smooth tapered jaws

**RS-6430**

CASTROVIEJO
round handle with lock
carbide jaws
length 4 1/2"



RS-6437 Shown

RS-6437

Micro Needle Holder
0.6mm tips
length 5 1/8"
curved jaws with lock



RS-6445 Shown

RS-6445

Micro Needle Holder
0.6mm tips
length 7 1/8"
curved jaws with lock



RS-7800

Intracardiac
extra delicate
with carbide jaws

length 5"

RS-7804

length 7"




Note: tungsten carbide inserts have 4,000 teeth per square inch


RS-7820

BABY DERF
large finger rings
length 4^{1/4}"



1-800-424-2984
www.roboz.com





RS-7822 Shown

RS-7822
RS-7824


DERF
length 4³/₄"
standard jaws
carbide jaws



RS-7830 Shown

RS-7830
RS-7832
RS-7842

WEBSTER
smooth jaws, delicate
length 4³/₄"
length 5"
serrated carbide jaws
length 5"





- RS-7840** HALSEY
length 5"
smooth jaws
- RS-7841** serrated jaws
carbide
- RS-7841L** left handed version



RS-7840 Shown

- RS-7850** HEGAR-
BAUMGARTNER
length 5 1/4"



RS-7860

RS-7862

CRILE-WOOD

length 6"

length 7"



RS-7860 Shown

RS-7880

RS-7882

OLSEN-HEGAR
extra delicate with
combination scissor
length 4^{3/4}"

serrated jaws

extra delicate

smooth



RS-7880 Shown



OLSEN-HEGAR
needle holder scissors
combination cross
serration tip

RS-7884

length 5 1/2"

RS-7886

length 6 1/2"



RS-7884 Shown

OLSEN-HEGAR
needle holder scissors
combination
with carbide jaws
length 4 3/4" TC


RS-7892

RS-7894

length 5 1/2"


RS-7896

length 6 1/2"




RS-7894 Shown

Note: tungsten carbide inserts have 2,600 teeth per square inch



RS-7900 METZENBAUM
length 7"
straight




RS-7910 MAYO-HEGAR
length 5"
RS-7912 length 6"
RS-7912L left handed version
RS-7920 with carbide jaws
length 5"
RS-7922 with carbide jaws
length 6"




RS-7910 Shown


Note: tungsten carbide inserts have 2,600 teeth per square inch




RS-7960
BROWN
convex jaw
with fine serrations
for easy suture tying
length 5¹/₄"



RS-7960 Shown



**Reflex Wound Clip Products**

- RS-9260** Reflex 9 Wound Clip
Applier for Reflex 9mm
clips only
- RS-9262** Reflex 9 Wound Clip
9mm for use with
RS-9260
box of 100
Stainless Steel
- RS-9265** Reflex 9 Wound Clip
9mm for use with
RS-9260
box of 1,000
Stainless Steel
- RS-9250** Reflex 7 Clip Applier
for Reflex 7mm clips
only
- RS-9258** Reflex 7 Clip 7mm
for use with RS-9250
box of 100
Stainless Steel
- RS-9255** Reflex 7 Clip 7mm
for use with RS-9250
box of 1,000
Stainless Steel



- RS-9263** Reflex Wound Clip
removing forceps
for 9mm, 7mm and
other brands
length 4"

1-800-424-2984
www.roboz.com



- MICHEL
wound clips
packaged 1,000 pcs.
per box
- RS-9270** 7.5mm
- RS-9272** 11mm



- RS-9290** HEGENBARTH
clip applying forceps
cross action
length 5"



- RS-9292** MICHEL
clip applying forceps
length 5"



- RS-9294** MICHEL
clip applying and
removing forceps
length 5"



COMPARISON OF MILLIMETERS TO INCHES

MM	INCHES	MM	INCHES
0.33	.013	6	.236
0.5	.020	6.33	.249
0.67	.026	6.5	.256
1	.039	6.67	.263
1.33	.052	7	.275
1.5	.059	7.33	.288
1.67	.066	7.5	.295
2	.079	7.67	.302
2.33	.092	8	.315
2.5	.098	8.33	.328
2.67	.105	8.5	.335
3	.118	8.67	.341
3.33	.131	9	.354
3.5	.138	9.33	.367
3.67	.144	9.5	.374
4	.157	9.67	.380
4.33	.170	10	.394
4.5	.177	10.33	.407
4.67	.184	10.67	.420
5	.197	11	.433
5.33	.210	11.33	.446
5.5	.216	11.67	.459
5.67	.223	12	.472

COMPARISON OF INCHES TO MILLIMETERS

INCHES	MM	INCHES	MM
1/32	0.79	17/32	13.49
1/16	1.59	9/16	14.29
3/32	2.38	19/32	15.08
1/8	3.18	5/8	15.88
5/32	3.97	21/32	16.67
3/16	4.76	11/16	17.46
7/32	5.56	23/32	18.26
1/4	6.35	3/4	19.05
9/32	7.14	25/32	19.84
5/16	7.94	13/16	20.64
11/32	8.73	27/32	21.43
3/8	9.53	7/8	22.22
13/32	10.32	29/32	23.02
7/16	11.11	15/16	23.81
15/32	11.91	31/32	24.61
1/2	12.7	1"	25.4

COMPARISON OF INCHES TO CENTIMETERS

INCHES	CM
1	2.54
2	5.08
3	7.62
4	10.16
5	12.7
6	15.24
7	17.78
8	20.32
9	22.86
10	25.40
11	27.94
12	30.48
13	33.02
14	35.56
15	38.10

COMPARISON OF CENTIMETERS TO INCHES

CM	INCHES
1	.3937
2	.7874
3	1.1811
4	1.5748
5	1.9685
6	2.3622
7	2.7559
8	3.1496
9	3.5433
10	3.937
20	7.8740
30	11.811
40	15.748

